Tucker-Davis Technologies TTank format:

Data collected or used by TDT software is stored in *tanks* that are logically divided into *blocks*. These blocks contain all the data from a single recording session such as an experimental run. Within each block different stores can record different types of events at different rates. This document is intended to provide the user with the details necessary to directly access data stored in these tanks.

The disk format for a TDT tank consists of a single folder (the tank directory) containing per-block subfolders. Hidden *desktop.ini* files are used to set the icon on all of these folders and as a means for TDT software to recognize tank folders. Within each block folder are four additional files each named according to "TANKNAME_BLOCKNAME.XXX" where XXX is one of the following extensions: TBK, TDX, TSQ, TEV.

The low level details of the TBK and TDX files are beyond the scope of this document. Briefly, the TBK contains block-level headers and offsets to the beginning of each second of data collection in the TSQ file. The TDX file contains epoch indexing information for faster access or filtering when dealing with epoch data. The remaining two files (the TSQ and TEV) contain the actual event headers and data meaning that usable TBK and TDX files can be reconstructed from these files with tools provided by TDT technical support if the TBK or TDX becomes lost or corrupted. The TSQ holds the header information for all the individual events whose data are in turn stored in the TEV.

The format of the TSQ file is simply a list of headers for each data segment or snippet stored in the TEV. Fields in the headers are little-endian, aligned at 32-bit boundaries (except for Sort Code), and are interpreted as follows:

Field	Туре	# Bytes	Comment
Size	long	4	Size of data, with header, in longs.
Type	long	4	Snip, stream, scalar, etc.
Code	long	4	4-character event name (cast as long)
Channel	unsigned short	2	
Sort Code	unsigned short	2	Sort code for snip data.
Time Stamp	double	8	Event start time in seconds.
Event Offset or Strobe	int64 or double	8	Offset in TEV file or raw value.
Data Format	long	4	float, long, etc.
Frequency	float	4	Sampling frequency

This results in each header being 10 words, or 40 bytes, long. This constant size allows for faster random access within the TSQ files. Data in the TEV file can be variable length, for example if different types of data were being recorded into the same block. Two fields also require additional explanation. The *type* field specifies the type of data this header contains or refers to (snippets, streaming, etc.) and the *data format* specifies the format of the data (float, int, etc.). The possible values are shown in the following 2 tables:

Type	Value in Hexadecimal
EVTYPE_UNKNOWN	0x00000000
EVTYPE_STRON	0x00000101
EVTYPE_STROFF	0x00000102
EVTYPE_SCALER	0x00000201
EVTYPE_STREAM	0x00008101
EVTYPE_SNIP	0x00008201
EVTYPE_MARK	0x00008801

Data Format	Value
DFORM_FLOAT	0
DFORM_LONG	1
DFORM_SHORT	2
DFORM_BYTE	3
DFORM_DOUBLE	4
DFORM_QWORD	5

The remaining pages contain example C code that was created access a tank in order to create a CSV file of the recorded data. Please feel free to copy this code and adjust it for your own uses.

```
// Demo program showing how to read tsg file and tev file
#include <stdio.h>
// tank file structure
// tbk file has block events information and on second time offsets
// to efficiently locate events if the data is queried by time .
// tsq file is a list of event headers, each 40 bytes long,
// ordered strictly by time .
//
// tev file contains event binary data
// tev and tsq files work together to get an event's data and
// attributes
//
// tdx file contains just information about epoc stores,
// is optionally generated after recording for fast retrieval
// of epoc information
/* names and paths used in the code snippet to locate the tank files */
#define TANK NAME "DEMOTANK2"
#define BLOCK NAME "Block-3"
#define EV NAME "LFPs"
#define CHANNEL 1
#define TANK PATH "C:\\TDT\\OpenEx\\Tanks\\" TANK NAME "\\" BLOCK NAME "\\"
/* TTank event header structure */
struct tsqEventHeader
{
    long
                   size;
                   type;    /* event type: snip, pdec, epoc etc */
code;    /* event name: must be 4 chars, cast as a long */
    long
                   type;
    unsigned short channel;  /* data acquisition channel */
unsigned short sortcode; /* sort code for snip data */
              timestamp; /* time offset when even occurred */
    double
    union
          int64 ev offset; /* data offset in the TEV file */
                  strobe; /* raw data value */
        double
    };
                              /* data format of event: byte, short, float
    long
                   format;
(usually), double */
    float
                  frequency; /* sampling frequency */
};
int main(int argc, char* argv[])
    FILE *tsq = NULL, *tev = NULL, *out = NULL;
    char bev[10000];
    tsqEventHeader bsq;
    tsq = fopen(TANK PATH TANK NAME " " BLOCK NAME ".tsq", "r+b");
    tev = fopen (TANK PATH TANK NAME " " BLOCK NAME ".tev", "r+b");
    out = fopen (TANK PATH TANK NAME " " BLOCK NAME ".csv", "w");
    if (tsq != NULL && tev != NULL && out != NULL)
        while(fread(&bsq, sizeof(bsq), 1, tsq) == 1)
```

```
if ( bsq.code == *(long*)EV_NAME && bsq.channel == CHANNEL)
{
    long nPts = bsq.size - 10;
    fsetpos(tev, &bsq.ev_offset);
    fread(bev, nPts <<2, 1, tev);
    for(long i = 0; i < nPts; i ++)
        fprintf(out, "%E,", ((float*) bev)[i]);
    fprintf(out, "\r\n");
}

if (tsq)
    fclose(tsq);
if (tev)
    fclose(tev);
if (out)
    fclose(out);

return 0;
}</pre>
```